

Date: 08/09/2025

Urgent Field Safety Notice

Various Clear-Therm™ Mini HMEF

For Attention of*: MDSO's, All clinical staff, Managers and Users of the above products, including those who may use these devices remotely.

Contact details of local representative (name, e-mail, telephone, address etc.)*

Chris Randall
Wokingham Site Quality Manager
Intersurgical Ltd.
Crane house
Molly Millars Lane
Wokingham
Berkshire
RG41 2RZ

Email: priority@intersurgical.co.uk
Tel. No: 0118 9656 300


Urgent Field Safety Notice (FSN)

Various Clear-Therm™ Mini HMEF

Risk addressed by FSN

1. Information on Affected Devices*																																															
1.	1. Device Type(s)*																																														
	Various Clear-Therm™ Mini HMEF																																														
1.	2. Commercial name(s)																																														
	<ul style="list-style-type: none">Clear-Therm™ Mini HMEF with luer portClear-Therm™ Mini HMEF with luer port and elbowEconomy, anaesthetic face mask, size 2, paediatric, 22F elbow and HMEF																																														
1.	3. Unique Device Identifier(s) (UDI-DI)																																														
	<ul style="list-style-type: none">1831000 - 50302670551591831197 - 50302670551971513001 - 5030267125470																																														
	4. Primary clinical purpose of device(s)*																																														
	Clear-Therm Mini HMEF is intended for reducing the risk of bacterial and viral contamination of patients, medical devices and equipment, whilst also reducing moisture and heat loss from the patient's respiratory gases within anaesthesia, critical and respiratory care breathing systems.																																														
1.	5. Device Model/Catalogue/part number(s)*																																														
	<ul style="list-style-type: none">183100018311971513001																																														
1.	6. Software version																																														
	N/A																																														
1.	7. Affected serial or lot number range																																														
	Any of these part numbers within the following Lot numbers. <ul style="list-style-type: none">1831000 <table><tr><td>32291126</td><td>32317717</td><td>32323746</td><td>32412583</td><td>32421306</td><td>32503642</td></tr><tr><td>32291130</td><td>32319201</td><td>32324614</td><td>32413321</td><td>32422521</td><td>32505952</td></tr><tr><td>32291139</td><td>32319998</td><td>32390307</td><td>32414436</td><td>32423219</td><td>32506067</td></tr><tr><td>32310682</td><td>32320554</td><td>32391058</td><td>32415119</td><td>32424642</td><td>32506824</td></tr><tr><td>32310945</td><td>32321330</td><td>32402792</td><td>32415685</td><td>32426677</td><td></td></tr><tr><td>32313487</td><td>32321698</td><td>32405756</td><td>32418358</td><td>32427421</td><td></td></tr><tr><td>32317296</td><td>32322273</td><td>32406802</td><td>32419265</td><td>32502751</td><td></td></tr></table>						32291126	32317717	32323746	32412583	32421306	32503642	32291130	32319201	32324614	32413321	32422521	32505952	32291139	32319998	32390307	32414436	32423219	32506067	32310682	32320554	32391058	32415119	32424642	32506824	32310945	32321330	32402792	32415685	32426677		32313487	32321698	32405756	32418358	32427421		32317296	32322273	32406802	32419265	32502751
32291126	32317717	32323746	32412583	32421306	32503642																																										
32291130	32319201	32324614	32413321	32422521	32505952																																										
32291139	32319998	32390307	32414436	32423219	32506067																																										
32310682	32320554	32391058	32415119	32424642	32506824																																										
32310945	32321330	32402792	32415685	32426677																																											
32313487	32321698	32405756	32418358	32427421																																											
32317296	32322273	32406802	32419265	32502751																																											

	<ul style="list-style-type: none">• 1831197<table><tr><td>32320639</td><td>32324900</td><td>32413019</td><td>32417583</td><td>32422536</td><td>32504025</td></tr><tr><td>32322207</td><td>32400440</td><td>32414189</td><td>32421768</td><td>32426712</td><td></td></tr></table>• 1513001<table><tr><td>32321802</td></tr><tr><td>32421512</td></tr></table>	32320639	32324900	32413019	32417583	32422536	32504025	32322207	32400440	32414189	32421768	32426712		32321802	32421512
32320639	32324900	32413019	32417583	32422536	32504025										
32322207	32400440	32414189	32421768	32426712											
32321802															
32421512															
1.	8. Associated devices														
	N/A.														

2. Reason for Field Safety Corrective Action (FSCA)*	
2.	1. Description of the product problem* <p>The two housings of some HMEF devices have been found to separate when a force is applied during handling of these products or movement and torsion due to repositioning of the patient. See separation of the housings shown below.</p> 
2.	2. Hazard giving rise to the FSCA* <p>The reported separation of the HMEF would cause gross leakage from the device, which in turn would result in leakage of gas from the breathing system. This could have the effect of the patient not receiving the prescribed mixture of anaesthetic gas and/or the prescribed ventilation, reducing the FiO₂ of inspired gases and resulting in the patient becoming hypoxic.</p>
2.	3. Probability of problem arising <p>Our investigation and inspection of potentially affected stock has estimated the probability of failure rate to be unlikely, which equals to 0.01% to 0.001% (1 in 10 000 to 1 in 100 000 products).</p>

2.	4. Predicted risk to patient/users The risks associated with the identified fault have been reviewed, and whilst the probability of occurrence is low, we believe it is essential to address the issue promptly to further reduce the risk of any potential patient harm or inconvenience to users.
2.	5. Further information to help characterise the problem N/A
2.	6. Background on Issue Following customer reports from the market and subsequent investigation including inspection of potentially affected stock and statistical analysis, we have determined that some products have been manufactured with inadequate ultrasonic weld of the two housings.
2.	7. Other information relevant to FSCA N/A
	3. Type of Action to mitigate the risk*
3.	1. Action To Be Taken by the User* <div style="display: flex; justify-content: space-between;"> <input checked="" type="checkbox"/> Identify Device <input checked="" type="checkbox"/> Quarantine Device <input type="checkbox"/> Return Device <input type="checkbox"/> Destroy Device </div> <div style="margin-top: 10px;"> <input type="checkbox"/> On-site device modification/inspection <input type="checkbox"/> Follow patient management recommendations <input type="checkbox"/> Take note of amendment/reinforcement of Instructions For Use (IFU) </div> <div style="margin-top: 10px;"> <input checked="" type="checkbox"/> Other <input type="checkbox"/> None </div> <p>Please distribute this Field Safety Notice to all potential users of the Clear-Therm™ Mini HMEF listed above, within your facility. This is for their awareness of the potential problem and to carry out the following actions.</p> <p>Please note: This is not a product removal due to the current limited availability and to allow continuity of supply of the Clear-Therm™ Mini HMEF devices.</p> <p>To ensure the safety of patients we recommend the following actions are taken in addition to those detailed in the IFU provided with the device.</p> <ol style="list-style-type: none"> 1. Identify any potentially affected products from the affected codes and lot numbers listed above. 2. Immediately before use unpack the device and carry out the check as described below. 3. All users must check the HMEF housing is securely welded by holding both connection tapers at each end and apply a downward breaking force by flexing the tapers in the direction shown by the arrows below.




4. Retain any affected sample(s) identified, and please report to us immediately.

Please complete and return the Reply Form provided to priority@intersurgical.co.uk, to confirm receipt of this notice and that the necessary actions are being taken.

Please continue to report to Intersurgical any adverse events involving this product.

3.	2. By when should the action be completed?	Immediately on receipt of this FSN, and awareness of this FSN should be ongoing until all potentially affected stock listed in this FSN has been used up.
3.	3. Particular considerations for: N/A Is follow-up of patients or review of patients' previous results recommended? Not applicable.	
3.	4. Is customer Reply Required? * (If yes, form attached specifying deadline for return)	Yes
3.	5. Action Being Taken by the Manufacturer <div style="display: flex; justify-content: space-between;"> <div> <input type="checkbox"/> Product Removal <input type="checkbox"/> Software upgrade <input checked="" type="checkbox"/> Other </div> <div> <input type="checkbox"/> On-site device modification/inspection <input type="checkbox"/> IFU or labelling change <input type="checkbox"/> None </div> </div>	

	Corrective actions have been implemented in the manufacturing process in April 2025 to eliminate this problem for all current and future supply.	
3	6. By when should the action be completed?	One month from receipt of the FSN
3.	7. Is the FSN required to be communicated to the patient /lay user?	No
3	8. If yes, has manufacturer provided additional information suitable for the patient/lay user in a patient/lay or non-professional user information letter/sheet?	
	N/A	
4. General Information*		
4.	1. FSN Type*	New – Advisory Notice
4.	2. For updated FSN, reference number and date of previous FSN	N/A
	3.	
4.	4. For Updated FSN, key new information as follows:	
	N/A	
4.	5. Further advice or information already expected in follow-up FSN? *	No
4	6. If follow-up FSN expected, what is the further advice expected to relate to:	
	N/A	
4	7. Anticipated timescale for follow-up FSN	N/A
4.	8. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	Intersurgical Ltd.
	b. Address	Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ
	c. Website address	https://www.intersurgical.com/
4.	9. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *	
4.	10. List of attachments/appendices:	Customer Reply Form
4.	11. Name/Signature	Ivan Seniut, Group Quality and Regulatory Affairs Director, Intersurgical
		

Transmission of this Field Safety Notice	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p>

Note: Fields indicated by * are considered necessary for all FSNs. Others are optional.

Field Safety Notice Customer Reply Form

1. Field Safety Notice (FSN) information																																																									
FSN Reference number*	491800																																																								
FSN Date*	08/09/2025																																																								
Product/ Device name*	<ul style="list-style-type: none"> • Clear-Therm™ Mini HMEF with luer port • Clear-Therm™ Mini HMEF with luer port and elbow • Economy, anaesthetic face mask, size 2, paediatric, 22F elbow and HMEF 																																																								
Product Code(s)	<ul style="list-style-type: none"> • 1831000 • 1831197 • 1513001 																																																								
Batch/Serial Number (s)	<p>Any of these part numbers within the following Lot numbers.</p> <ul style="list-style-type: none"> • 1831000 <table border="1"> <tbody> <tr><td>32291126</td><td>32317717</td><td>32323746</td><td>32412583</td><td>32421306</td><td>32503642</td></tr> <tr><td>32291130</td><td>32319201</td><td>32324614</td><td>32413321</td><td>32422521</td><td>32505952</td></tr> <tr><td>32291139</td><td>32319998</td><td>32390307</td><td>32414436</td><td>32423219</td><td>32506067</td></tr> <tr><td>32310682</td><td>32320554</td><td>32391058</td><td>32415119</td><td>32424642</td><td>32506824</td></tr> <tr><td>32310945</td><td>32321330</td><td>32402792</td><td>32415685</td><td>32426677</td><td></td></tr> <tr><td>32313487</td><td>32321698</td><td>32405756</td><td>32418358</td><td>32427421</td><td></td></tr> <tr><td>32317296</td><td>32322273</td><td>32406802</td><td>32419265</td><td>32502751</td><td></td></tr> </tbody> </table> <ul style="list-style-type: none"> • 1831197 <table border="1"> <tbody> <tr><td>32320639</td><td>32324900</td><td>32413019</td><td>32417583</td><td>32422536</td><td>32504025</td></tr> <tr><td>32322207</td><td>32400440</td><td>32414189</td><td>32421768</td><td>32426712</td><td></td></tr> </tbody> </table> <ul style="list-style-type: none"> • 1513001 <table border="1"> <tbody> <tr><td>32321802</td></tr> <tr><td>32421512</td></tr> </tbody> </table>	32291126	32317717	32323746	32412583	32421306	32503642	32291130	32319201	32324614	32413321	32422521	32505952	32291139	32319998	32390307	32414436	32423219	32506067	32310682	32320554	32391058	32415119	32424642	32506824	32310945	32321330	32402792	32415685	32426677		32313487	32321698	32405756	32418358	32427421		32317296	32322273	32406802	32419265	32502751		32320639	32324900	32413019	32417583	32422536	32504025	32322207	32400440	32414189	32421768	32426712		32321802	32421512
32291126	32317717	32323746	32412583	32421306	32503642																																																				
32291130	32319201	32324614	32413321	32422521	32505952																																																				
32291139	32319998	32390307	32414436	32423219	32506067																																																				
32310682	32320554	32391058	32415119	32424642	32506824																																																				
32310945	32321330	32402792	32415685	32426677																																																					
32313487	32321698	32405756	32418358	32427421																																																					
32317296	32322273	32406802	32419265	32502751																																																					
32320639	32324900	32413019	32417583	32422536	32504025																																																				
32322207	32400440	32414189	32421768	32426712																																																					
32321802																																																									
32421512																																																									

2. Customer Details	
Account Number	
Healthcare Organisation Name*	
Organisation Address*	
Department/Unit	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

3. Customer action undertaken on behalf of Healthcare Organisation				
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understood its content.	Customer to complete or enter N/A		
<input type="checkbox"/>	The information and required actions have been brought to the attention of all relevant users and executed.	Customer to complete or enter N/A		
<input type="checkbox"/>	I do not have any affected devices.	Customer to complete or enter N/A		
<input type="checkbox"/>	We have the following potentially affected stock we wish to return for credit/replacement. (Please enter the quantity for each Code and Lot number).	Code:	Lot:	Qty:
		Code:	Lot:	Qty:
		Code:	Lot:	Qty:
		Code:	Lot:	Qty:
		Code:	Lot:	Qty:
<input type="checkbox"/>	Any Other comments:			
Print Name*		Customer print name here		
Signature*		Customer sign here		
Date*				
4. Return acknowledgement to sender				
Email		priority@intersurgical.co.uk		
Customer Helpline		N/A		
Postal Address		Intersurgical Ltd., Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ		
Web Portal		N/A		
Deadline for returning the customer reply form*		10/10/2025		

Mandatory fields are marked with *

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN.

Your organisation's reply is the evidence that the manufacturer, and the MHRA, needs to monitor the progress of the corrective actions to ensure patient safety.

Without your reply the manufacturer can't know if their important message has been received and the MHRA may need to issue a safety communication.