

Quick-Stop by Bray Silver Nitrate Applicator

Title

Quick-Stop 75 % Silver Nitrate Applicator by Bray – Flexible Single-Use Nail, Paw & Minor Wound Care for Pets – Veterinary-Grade Tool You Can Use At Home for Grooming & Small Injuries

Key Features / Short Description

- **Precise Minor Wound Care:** Single-use 75 % silver nitrate applicator designed for pet nail and small surface injuries, ideal when a nail is overly-trimmed and the quick is exposed or on a bleeding wound.
- **Flexible & Easy to Use:** Bendable 15 cm plastic handle provides accurate application to the pet's nail or minor scrape.
- **Hygienic Single-Use Design:** Each applicator is packaged for convenience, with option for multiple to used if required.
- **Trusted Veterinary Quality:** Natural antimicrobial action supports clean, healthy recovery for pets during grooming or first-aid care.
- **Portable & Ready-to-Go:** Slim, lightweight applicator fits first-aid kits or grooming bags for on-the-spot pet care.

Long Description

The **Quick-Stop Silver Nitrate Applicator by Bray** delivers professional-grade first-aid care for small pet injuries such as an **over-clipped nail** or other **grazes or cuts**. Whether you're a groomer, vet, or attentive **pet parent**, this 75 % **silver nitrate** applicator provides precise, single-use care when a **pet gets injured**.

Using the applicator is simple. Gently bend the 15 cm plastic handle for better control, moisten the tip with a drop of water or allow natural moisture from the bleeding wound to trigger the **silver nitrate**, and carefully touch the **injury**—no need to **apply pressure**. If needed, **remain calm**, hold a clean cloth to the injury to support natural clotting, and dispose of the applicator after use. The **single-use** format ensures hygiene and convenience for both **pet owners**, grooming professionals, and veterinarians in the practise or during surgery.

Dogs' nails contain a sensitive quick rich in **blood vessels** and nerves. A misjudged clip or torn **dew claw** can expose the **nail bed**, causing discomfort and natural bleeding. A **styptic applicator** helps protect the **dog's foot** and encourages a protective **blood clot** to form so the area is quickly managed. Compared with ordinary **styptic powder**, this applicator offers accurate, mess-free

application directly to the **pet's injury**, and with antimicrobial/antibacterial properties helping the wound stay clean while nature does the rest.

Routine nail care remains important. Keep your **dog's nails** and **black nails** neatly **trimmed** to reduce the risk of a painful **broken nail**. During regular **nail clipping**, trim a little at a time to avoid the sensitive quick and maintain a healthy **nail bed**. Even with care, the **bleeding quick** can be reached; having this **first aid** essential on hand means you're prepared.

Beyond nail care, the Quick-Stop Silver Nitrate Applicator can also help manage other small pet injuries. Minor scrapes, tiny surface cuts, or nicks from grooming tools on areas such as paws, ears, or tails benefit from the applicator's precise application and natural antimicrobial action, providing quick, convenient support for everyday pet mishaps.

Made in the UK to exacting standards, each applicator is individually packaged for sterility and easy transport. Store in a cool, dry place until needed. Intended for animal use only and not licensed for humans, it's a trusted choice for veterinarians, groomers, and **pet parents** who want clean, effective care for minor wounds and nail-related incidents.

Safety Notes

- For animal use only – not for humans.
- Contains corrosive silver nitrate & potassium nitrate; avoid skin/eye contact and protect surfaces (may stain).
- Not for large or surgical wounds.
- Store in original packaging away from light and moisture.



