


The ideal thread for the uterine suture



Surgicryl® Monofast Polyglecaprone 25
about 20 days of tissue support
complete resorption after 90-120 days

- ☑ absorbable synthetic
- ☑ monofilament
- ☑ resorption by hydrolysis

☑ Directions:

soft tissues, uterus, bladder, intestines, etc.

SURGICRYL® MONOFAST has the following advantages:

- ☑ faster resorption of SURGICRYL® MONOFAST with better healing of uterine tissue
- ☑ good tensile strength
- ☑ Excellent maneuverability due to the high flexibility of the wire
- ☑ No saw effect due to its flexible monofilament structure
- ☑ Excellent stability with a new preservative liquid



The faster resorption of SURGICRYL® MONOFAST has a direct influence on fertility, a Factor of prime importance in the profitability of cattle production.
SURGICRYL® MONOFAST is an important advance for modern veterinary medicine.

ON WHEELS

Metric	USP	Length	Code
4	1	150 cm	16400015
5	2	150 cm	16500015
5	2	250 cm	16500025

